

Safety Trays



Safety Trays and Liners

Unlike disposable paper bench protectors, Scie-Plas safety trays and tray liners are reusable. This reduces running costs and drastically reduces contaminated waste volume. The rigid PVC base features specially designed stabilising edges and rounded corners for easy cleaning. Additionally, the liner's non-porous surface allows valuable samples to be retrieved - something which is clearly not possible with absorbent paper protectors.

- Saves on running cost
- Contains liquid spillages
- Provides clearly defined work area
- Easily cleaned
- APET - environmentally friendly when incinerated
- Suitable for most hazardous spillages
- Anti-static and non-porous

ORDERING INFORMATION

	Tray Type	Ext Dimensions	Int Dimensions	Base only	APET Liners pk/25
 General Purpose		46 x 26cm	40.5 x 20cm	RPP-TW4626	RPP-TL4626
		54 x 34cm	46 x 26cm	RPP-TW5434	RPP-TL5434
		57 x 54cm	45.5 x 43cm	RPP-TW5754	RPP-TL5754
		68 x 54cm	56.5 x 42.5cm	RPP-TW6854	RPP-TL6854
		70 x 46cm	57 x 35cm	RPP-TW7046	RPP-TL7046
		113 x 54cm	100 x 42cm	RPP-TW11354	RPP-TL11354
		 Radiation Hazard		46 x 26cm	40.5 x 20cm
54 x 34cm	46 x 26cm			RPP-TY5434	RPP-TL5434
57 x 54cm	45.5 x 43cm			RPP-TY5754	RPP-TL5754
68 x 54cm	56.5 x 42.5cm			RPP-TY6854	RPP-TL6854
70 x 46cm	57 x 35cm			RPP-TY7046	RPP-TL7046
113 x 54cm	100 x 42cm			RPP-TY11354	RPP-TL11354
 Biohazard				46 x 26cm	40.5 x 20cm
		54 x 34cm	46 x 26cm	RPP-TO5434	RPP-TL5434
		57 x 54cm	45.5 x 43cm	RPP-TO5754	RPP-TL5754
		68 x 54cm	56.5 x 42.5cm	RPP-TO6854	RPP-TL6854
		70 x 46cm	57 x 35cm	RPP-TO7046	RPP-TL7046
		113 x 54cm	100 x 42cm	RPP-TO11354	RPP-TL11354